




## ***ACHERON INSTRUMENTS***

Superior quality and attention to detail make our **Premium** a long-term value for your OR.

- Exceptional warranties ensure your instruments continue to perform at a high level
- Crafted using 400 series grade stainless steel for increased durability
- Stellite inserts in our scissors create a sharper, smoother cutting action
- Tungsten Carbide jaws are designed to sustain their straight edge and needle-holding power



***Premium & Economy Grade  
Instruments***


# Acheron Instruments

Our Premium grade instruments represent our finest surgical instrumentation. The jaws on instruments are precisely ground and shaped for enhanced and accurate performance, while the satin finish is designed to reduce glare during surgical procedures.

## Needle Holders

Tungsten Carbide Inserts in our needle holders are guaranteed for **3 years** from date of purchase.

- Tungsten Carbide jaws are designed to sustain their straight edge and needle-holding power for years
- Milled inserts provide uniform serrations to enhance stability when grasping a needle
- Tips are uniformly smooth and carefully rounded to avoid catching sutures
- Jaw serrations are precisely meshed for a firmer grip
- Rachets are precisely undercut and polished for secure, effortless locking and unlocking
- Shanks are slim and lightweight for easy, flexible use



## Scissors

Our premium scissor inserts are guaranteed to stay sharp for **3 years** from date of purchase under normal use.

- TC inserts in our Scissors create a sharper, smoother cutting action.
- Tips are uniformly rounded on all blunt scissors and are designed for precise distal cutting
- Cutting edges with Tungsten Carbide inserts stay sharp longer
- Blades are carefully contoured to provide precise scissoring action



## Forceps

TC Inserts in our forceps are guaranteed for **3 years** from date of purchase.

- Tungsten Carbide jaws allow a firm hold with light pressure
- Tips double as both tooth and smooth because they will hold securely without penetrating the tissue
- Jaw serrations are precisely meshed to provide firmer grip



# Acheron Instruments: Economy Grade

Our instruments are a culmination of more than 30 years of experience in the design and manufacture of surgical instruments. They are made of highly-tempered stainless steel, which combines rust and corrosion resistance with rugged, long-lasting strength. All **Economy Grade** instruments are unconditionally **guaranteed for 6 months**.

## Hemostats

Our **Economy Grade** hemostats are guaranteed for **6 months** from date of purchase under normal use.

- Flat screws provide a smooth, clean profile to prevent tissue abrasion
- Tips are uniformly smooth and carefully rounded to prevent snagging
- Box lock on inner surfaces allows for a smooth and effortless action, while cold assembly followed by tempering adds extra strength and durability
- Shanks are slim and lightweight for easy, flexible use



## Scissors

Our **Economy Grade** scissors are guaranteed to stay sharp for **6 months** from date of purchase under normal use or we will sharpen them for free.

- Flat screw provides a smooth, clean profile to prevent tissue abrasion
- Finish is a unique combination of bright and matte to help distinguish from other instruments



Visit : [www.Acheron-Instruments.com](http://www.Acheron-Instruments.com) for more information.

## Warranty Details

Acheron Instruments offers specified Warranty Programs as outlined below:

Premium branded surgical instruments are warrantied for workmanship and materials for five (5) years. Should anything happen to disrupt the function of these instruments, we will repair them at no charge to you.

All Tungsten Carbide inserts on Premium Plus Needle Holders & Forceps are guaranteed for five (5) years of normal use. All Stellite inserts in Premium Plus Scissors are guaranteed for ten (10) years of normal use.

All Acheron instruments are warranted to be free from defects in workmanship and materials or will be replaced or repaired at no charge, at Acheron Instruments discretion. Warranty assumes product maintenance, if required, has been preformed at Acheron Instruments authorized repair center. Alteration, repair or modification of any product performed by persons not authorized by Acheron Instruments will result in immediate loss of warranty.

# Acheron Instruments

**Production Unit:** Industrial Estate Shahab Pura, Ugoki Road, Near Saga Chowk Sialkot, Pakistan.

**Head Office:** Building # 349 Block A, Jubilee Town Canal Road Lahore Pakistan | +923334301417 | [info@Acheron-Instruments.com](mailto:info@Acheron-Instruments.com) |

[www.Acheron-Instruments.com](http://www.Acheron-Instruments.com)

©2020 Acheron Instruments.. All rights reserved.